

**Active Solutions for Health** 

# Advances in peptide synthesis with green chemistry and cutting-edge technologies

Oleg Babii, Eugenie Fournier, Arjun Vijeta, Mohamad-Jamal Wawi, Bernd Henkel, Hélène Adihou

## Partner with EUROAPI

- Extensive experience with highly modified oligonucleotides & peptides, challenging chemistries and sequences
- Time- and cost-efficient transition from small to commercial scale production
- Investment in capability, process optimization & innovation
- Comprehensive support for the preparation of regulatory CMC package
- Strong customer interaction and transparent project management

# Innovation for SPPS

### Several existing challenges in SPPS:

- Green chemistry
- DMF restriction, hazardous to the environment & human health
- Large excess of reagents needed; sometimes a coupling needs to be repeated
- Substantial solvent consumption, especially for the washing steps
- Process development time: complex processes with multiple steps

## Greener solvents

Development of new green solvent mixtures to replace DMF in peptide synthesis: <u>Ternary mixture</u> (CMR free, class 3\*)

→ Successful synthesis of Bivaluridin (20 mer), Thymosin alpha (29 mer) and a GLP1 agonist (39 mer)


Origin of solvent mixture	Solvents	Tested ratio	Thymosin a purity (%)	Yield (%)
Standard	DMF	-	28	58
Pawlas et al. Schönleber, Pedersen et al.	DMSO/EtOAc	3:7	25	59
Pawlas et al.	NBP/EtOAc	1:3	51	69
Papini, Pacini, Rovero et al.	DMSO/BuOAc	3:7	23	58
EUROAPI	Ternary mixture	-	53	80

→ Enhance purity and yield of peptides while using class 3 solvents (low toxic potential to man) to promote sustainable development processes

## Reduction of solvent consumption

Current SPPS process uses batch washes which require large volume of solvent and long operation time.


Simulated flow through approach for wash steps is proposed.




- 60% solvent consumption reduction due to more efficient wash process.
- 20% operation time reduction due to faster wash process in flow mode.
- Improved crude purity due to shorter holding time between coupling steps.
- Low regulatory impact. Drug Muster File does not describe wash process because during washes not chemical reaction occurs.
- Minimal CAPEX investment would be required.
- Timeline for lab scale POC 6 weeks.
   Successful POC: Thymosin alpha (29 mer), GLP1 agonist (39 mer)
- Overall cost saving due to increased process intensity and reduce solvent use.

# Real time monitoring and Rotating bed reactor technology for SPPS

# New experimental set-up for online monitored SPPS





## Aib-enkephalin

Classical short model peptide
Challenging synthesis: steric
hindrance due to Aib
residues

2400		A PORT		10.1	492	- 450		600 -	•							
2200		September 1	1 4			- 425 em ox	Fmoc	500						•		•
1800 -	13.	Salah Salah Salah	-	and the same of th		400 A	Intensity @Fmoc	400			-					
1600					400-00	- 375 <u><b>å</b></u>										
1400						- 350		300				port	AP.	a dead		
1200				1		- 325		200	9							
1000	20			60	80	1 300		100	5	1			20	25	30	
		Time	(minute)				(d)				III	me (mir	nutej			
(c)																
		Deprote	ection kin	etics		700		ı.	_		Depre	otection	n wash	ning		
(c)	akur Çanê	. p.;	ection kind	etics	. 129	700		8			Depre	otection	n wash	ning	•	
6500 <del>-</del> 6000 -	g patrof, mi	. p.;	ection kind	etics	ne projet			8-		•	Depre	otection	n wash	ning	٠	
6500 <del>-</del> 6000 -	gareriae	. p.;	ection kin	etics	and the second	- 600 - 500			-	•	Depr	otection	n wash	ning	•	
6500 <del>-</del> 6000 -	garateria	. p.;	ection kind	etics	الاستان الاستان	- 600 - 500 - 400 @				•	Depr	otection	n wash	ning	•	
erdine 6000 -	godere energy	. p.;	المالية المالية المالية	etics	anica)	- 600 - 500 - 400 @		6-			Depre	otection	n wash	ning	•	
6500 <del>-</del> 6000 -	gentre see	. p.;	المالية المالية المالية	etics	and the second	- 600 - 500 - 400 @ - 300	Piperidine concentration(%)	6-			Depro	otection	n wash	ning		0.

(a) coupling kinetics of Fmoc-Phe-OH(b) coupling washing of Fmoc-Phe-OH(c) deprotection kinetics of Fmoc-Aib-OH (4th amino acid)(d) deprotection washing of Fmoc-Aib-OH (4th amino acid)

Investigating the responses of Raman spectrometry and Refractive index for all steps:

Raman spectroscopy monitoring of reaction and washing steps during Aib-enkephalin synthesis

POC with synthesis of Gonadorelin (approved API):

Crude Gonadorelin purity	70 %
Purified Gonadorelin purity	91 %
Overall yield	60 %

### Major results:

Instant flagging of deviation
Quick decision making
Enhanced efficiency
Improved quality
Improved PMI

### Conclusion

PAT is a promising solution for accelerated peptide synthesis development

SpinChem is an attractive alternative technique to perform peptide synthesis

PAT is also applicable for standard peptide synthesis and oligonucleotides

		Pu	rity
	PMI kg/kg	Aib-enkephalin H-Tyr-Aib-Aib-Phe-Leu-NH2	DesAib-impurity H-Tyr-Aib-Phe-Leu-NH2
Optimized protocol SpinChem reactor	2541	78 %	15 %
Standard protocol CSBio synthesizer	3795	20 %	77 %